

## Cleaning and Disinfection

90cm wide Freestanding/built-under units

### G 7882CD Disinfector



### Performance specifications

AN sets: 2/E 435 TA  
DIN mesh trays: 4/E 327 plus 4 E 142

### Control system

MTN - Multitronic Novo Plus  
Electronic controls with 10 programmes  
Various process parameters can be modified  
Display of programme sequence and faults  
Temperature and countdown display

### Features

Stainless-steel cabinet with 2 levels  
Spray arms/Stainless steel: 2  
Docking system for upper baskets/mobile injector units  
Quadruple filtration system  
UP Qmax circulation pump 400 l/min  
AP Drain pump  
DK Steam condenser/aerosol injection  
ET Electric door lock  
PASI Programme interruption backup  
SU Buzzer (acoustic signal)  
SST Interface for process documentation  
TA Drying unit/radial fan  
WES Water softener

### Heating

Electric 400V 9kW conv to 230V 6kW

### Dispensing technology

1 door dispenser for powder detergent  
1 door dispenser for liquid products  
1 DOS 10/30 dispenser pump for liquid acidic agent  
1 dispenser pump 60/30 for liquid detergent  
DOS module connection: DOS G 60, dispenser module for liquid detergent

### Cleaning technology

Fresh water system  
Thermal disinfection at 93°C with 10 min. holding time  
Docking system for upper baskets and mobile injector units

### Design

Frontloading unit with bottom-hinged door  
Freestanding/built-under unit without lid  
Slanted fascia panel  
Outer casing:  
AE Stainless-steel  
H xW xD: 820 x 900 x 700 mm

### Test certificates

VDE, EMV/Interference suppression, DVGW,  
CE 0366, § 18 IfSG, WRAS